LB3013IMFRGBV2

IMPEDANCEMETRES GAMME STARTY

MODE D'EMPLOI

BODYFAT SCALES STARTY RANGE

INSTRUCTIONS

Caractéristiques du produit / Product Information / Característica del producto / Caratteristica del prodotto / Productkenmerk / Produkteigenschaft / Característica do produto / Produktens egenskaper / Produktfunksjon / Tuotteen ominaisuudet / Produktkarakteristik			
Type d'équipement Equipment Type Tipo de equipo Tipo di apparecchiatura Type apparaatuur Gerätetyp Tipo de equipamento Typ av utrustning Utstyrstype Laitetyyppi Type udstyr	Équipement adaptatif non-FHSS non-FHSS adaptive equipment Equipo adaptativo no-FHSS Apparecchiatura adattiva non-FHSS Adaptieve apparatuur niet-FHSS Adaptives Gerät, nicht-FHSS Equipamento adaptativo não FHSS Icke-FHSS-anpassad utrustning Ikke-FHSS adaptivt utstyr Ei-FHSS-mukautuva laite Adaptivt udstyr uden FHSS	Gain d'antenne Antenna Gain Ganancia de antena Guadagno dell'antenna Antenneversterking Antennengewinn Ganho da antena Antenneforsterkning Antennförstärkning Antennivoitto Antenneforstærkning	-0,61 dBi
Technologies prises en charge Supported Technologies Tecnologías compatibles Tecnologie supportate Ondersteunde technologieën Unterstützte Technologien Tecnologias suportadas Tekniker som stöds Støttede teknologier Tuetut teknologiat Understøttede teknologier	Bluetooth Low Energy	Catégorie du récepteur Receiver Category Categoría del receptor Categoría del ricevitore Ontvangercategorie Empfängerkategorie Categoria do recetor Mottakerkategori Mottagarkategori Vastaanottimen luokka Modtagerkategori	3
Plage de fréquences de fonctionnement Operation Frequency Range Rango de frecuencias de funcionamiento Gamma di frequenze operative Werkfrequentiebereik Betriebsfrequenzbereich Faixa de frequências de funcionamento Driftsfrekvensområde Driftfrekvensområde Käyttöaajuusalue Driftsfrekvensområde	2.4G ISM band (2402MHz-2480MHz)	Alimentation Power Supply Alimentación Alimentazione Voeding Stromversorgung Alimentação Strømforsyning Strömforsörjning Sähkönsyöttö Strømforsyning	batterie battery batería tramite batteria batterij DC 4,5 V / batterie bateria batteri batteri akku batteri
Type de modulation Modulation Type Tipo de modulación Tipo di modulazione Modulatietype Modulationstyp Tipo de modulação Moduleringstype Modulationstyp Modulaatiotyyppi Modulationstype	GFSK	Capacité de géolocalisation Geo-location capability Capacidad de geolocalización Capacità di geolocalizzazione Geolocatiecapaciteit Geolokalisierungsfähigkeit Capacidade de geolocalização Geolokalisering Geolokaliseringskapacitet Paikannuskyky Geolokaliseringskapacitet	Non No No No Nee Nein não nei nej ei nej
Débit de données Data Rate Velocidad de datos Velocità di trasmissione dati Gegevensnelheid Datenrate Taxa de transferência de dados Datagjenomstrømning Datagenomstrømning Tiedonsiirtonopeus Datahastighed	1 Mbps	Nombre de canaux Number of Channels Número de canales Numero di canali Aantal kanalen Anzahl der Kanäle Número de canais Antall kanaler Antal kanaler Antal kanaler Antal kanaler	40 canaux 40 Channels 40 canales 40 canali 40 kanalen 40 Kanäle 40 canais 40 kanaler 40 kanaler 40 kanaler 40 kanaler
Largeur nominale de canal Nominal Channel Bandwidth Ancho nominal del canal Larghezza nominale del canale Nominale kanaalbreedte Nominale Kanalbreite Largura nominal do canal Nominell kanalbredd Nominell kanalbredd Kanavan nimellisleveys Nominel kanalbredd	1 MHz	Puissance RF maximale en sortie Maximum RF Output Power Potencia máxima de salida RF Potenza RF massima in uscita Maximale RF-uitgangsvermogen Maximale HF-Ausgangsleistung Potência máxima de saída de RF Maksimal RF-utgangseffekt Maximal RF-uteffekt Maximal RF-uteffekt Maksimal RF-udgangseffekt	-4,34 dBm (E.I.R.P)

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PRECAUTIONS FOR USE – SAFE USE

Safety precautions

People with electronic medical implants, such as pacemakers, should not use this product, as it emits a low-intensity electrical current that travels through the body and may interfere with the functioning of such devices.

Pregnant women should only use the weight measurement function. The other functions are not designed for use by pregnant women.

Warning

It is prohibited to:

- use the device on persons with pacemakers or other electronic medical implants. The electrical current transmitted by the product may interfere with these devices and cause serious consequences.
- use this device if an error is detected before or during measurement. Use under these conditions may cause injury or incorrect results.
- Interpreting the data without professional advice. The evaluation of measurements and the implementation of dietary or exercise programmes must be carried out by a specialist. Self-analysis may present health risks.
- Allowing a child to carry the product. They could drop it and injure themselves. Keep it out of the reach of young children.
- Heat the product, batteries or battery pack. This may cause deformation, fire or injury.
- Stand on the edge of the platform. Do not jump on the platform. This increases the risk of falling and serious injury.
- Place the product on a slippery surface, such as a wet floor. You may slip or fall and seriously injure yourself. Water may enter the device and damage its internal components.
- Do not step on the display screen or control panel.
- Do not insert your fingers into the openings of the device.
- Do not use rechargeable batteries. Do not use different types of batteries at the same time or use old batteries with new ones. In the case of batteries, battery fluid may leak or the batteries may overheat and disintegrate, causing damage to the device or injury to the user.

This device is designed for domestic use. It is not intended for professional use (hospitals, medical practices, gyms). Intensive or professional use will void the warranty.

Precautions for use

- Install the product on a flat, solid surface, away from vibrations to ensure reliable measurements.
- Avoid violent shocks or vibrations that could damage the product.

- Do not step onto the platform with wet feet or body. Dry yourself before use to avoid falls and ensure accurate measurements.
- Do not use the device in extreme environments: direct sunlight, proximity to heat sources, significant temperature variations or vibrations.
- Do not allow the device to come into contact with water.
- Never immerse the device in water or clean it with aggressive products (benzene, solvents). Use a soft, slightly damp cloth with a neutral detergent, then wipe with a dry cloth.
- Do not disassemble or modify the device.
- Ensure that the batteries are inserted with the correct polarity and in the correct direction.
- If the batteries are not inserted correctly, battery fluid may leak and damage the product. It is advisable to remove the batteries before storing the device if you are not going to use it for a long period of time.
- Keep the device away from sources of electromagnetic waves (e.g. microwaves) to avoid interference.
- This device is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or persons with no experience or knowledge, unless they have been given supervision or prior instruction concerning the use of the device by a person responsible for their safety. Children should be supervised to ensure that they do not play with the device.
- Remove your socks or stockings before measurement and ensure that your feet are clean.
- Do not move during measurement to avoid errors. Keep your arms and legs apart during measurement.

Safety instructions for the use, storage and disposal of batteries:

- Keep batteries away from children
- Insert the batteries respecting indicated polarities
- Do not mix batteries with different compositions
- Do not mix batteries from different brands
- Do not mix new and used batteries
- In case of prolonged non-use, remove the batteries to prevent them from leaking
- Store and use the batteries in a temperature range between 10°C and 40°C
- Misuse can lead to serious risks (explosion, excessive heat, fire, smoke) or a reduction in their lifespan
- Batteries are not household waste and should not be thrown away with household waste: dispose of them in specific collection bins for this purpose or contact your local administration to find out how to throw them away without any risk to the environment.

Little Balance disclaims all responsibility in the event of non-compliant use, particularly in case of voltage higher than that indicated. Misuse can result in injury or property damage.

SPECIFICATIONS

- Units of measurement: kg/lb (changeable on the back of the device or via app)
- Weight limit: 5-180 kg
- Power: 3 AAA batteries
- Number of accounts: no limit – max. 24 users per account
- Operating temperature: 10-40°C
- Operating humidity: 20-90%

Following operating alerts may appear:

- Lo: low batteries – please replace
- Err: overload
- C: error, please step off the scale and step back on

You will find the following in the package:

- Bodyfat scale
- Instructions for use
- 3 AAA batteries

Thank you for placing your trust in us!

Our bodyfat scale will help you easily measure your key body parameters:

1. Weight
2. BMI (body mass index)
3. Body fat
4. Muscle mass
5. Lean body mass
6. Subcutaneous fat
7. Visceral fat
8. Body water
9. Skeletal muscle
10. Bone mass
11. Protein
12. BMR (basal metabolic rate)
13. Body age
14. Fat mass
15. Water weight
16. Protein mass
17. Ideal body weight
18. Obesity level
19. Body type

Scan and upload

iOS 8.0 ou +

Android 6.0 ou +

INSTRUCTIONS FOR USE

The Fitdays+ app is extremely intuitive and offers a choice of multiple languages (Account menu), which is why we have endeavoured to write short and simple instructions.

1. Download the Fitdays+ app from the Apple Store or Google Play; the corresponding icon will appear on your phone

2. Remove the battery insulation
3. Place the device on a hard, flat, dry surface (no carpet, rugs, etc.)
4. Ensure that Bluetooth is enabled on your phone.
5. Register with your email address or social media account.
6. Enter your personal details: please note that incorrect information will lead to inaccurate analysis results.
7. Activate the impedance meter by gently pressing down with your foot (the screen will light up).
8. Connect your phone via Bluetooth by selecting the following options:
 - Device
 - Add a device
 - My scale (symbolised by a personal scale)
 - you can rename your device if you wish
9. Before weighing yourself, check that the scale screen is on, or activate it by gently pressing down with your foot
10. Also check that your phone is turned on and connected, which you can check on the home page (otherwise, the scale will only display your weight on its screen)
11. Wait until the screen displays 0
12. Step onto the scale barefoot and with dry feet (if you step on with socks, the app will only be able to display your weight and calculate your BMI)
13. The screen will flash during the analysis and calculation process; remain still until your phone has received all the information from the scale
14. All analyses appear on your screen and can be scrolled through
15. Tapping on a data display opens a positioning relative to a standard value (please note: this is only an indicative positioning and does not replace medical advice)
16. Tap the 'trend' menu to see the evolution of each parameter over weekly periods, which you can scroll through to view over a longer period
17. When adding another user, the device performs automatic recognition and reminds you to check that the application is set to the correct profile
18. Baby mode:
 - On the home page, scroll down to baby mode.
 - Step onto the device alone and it will display your weight.
 - Then pick up your baby and the device will display your baby's weight.

TROUBLESHOOTING

😞 bluetooth won't connect :

- check that Bluetooth is enabled on both devices
- check your iOS version, which must be at least 8.0, or 6.0 on Android
- click on 'Device' on the home page
- then click on 'Add a device' in the 3 dots in top right-hand corner

😞 the analyses are not being performed :

- you must be barefoot to perform the analyses
- your feet (and the device) must be dry to perform the analyses
- wait until the numbers have finished flashing; this may take 10 seconds
- check the Bluetooth connection
- check that you are connected to the correct device

😞 the weight displayed is inconsistent :

- check that you are standing on a flat, dry and hard surface
- check the underside of the device to ensure that nothing is blocking any of the sensors
- remove the batteries and restart

😞 the analyses are incorrect :

- check that you have entered the correct body data
- check that you are logged in as the correct user

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Little Balance SAS - RCS de Versailles n° 798 750 246