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LB8662 Bodyfat scale - Operating instructions

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Thank you for placing your trust in Little Balance.

Although our equipment is very intuitive to use, we have endeavoured to provide clear and comprehensive instructions.

IMPORTANT: Before using the appliance, please read the instructions in this manual carefully. Failure to comply with the instructions or this manual may result in personal injury or damage to the appliance.

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CONSUMER SERVICE

You can contact Littlebalance at balancy@littlebalance.fr

The development team can be reached via the app: Account / comments / scales

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A screenshot of the Fitdays app registration screen. At the top, it shows the time '08:42' and signal strength. Below that, there are two buttons: 'Register' and 'Feedback'. The main form has two sections: 'Email' with a text input field and a 'Please enter your Email' prompt, and 'Password' with a text input field and a 'Please input password' prompt. Below the password field, there are two bullet points: 'password should contain 8-16 characters' and 'Include at least two kinds among numbers, uppercase and lowercase letters, symbols'. At the bottom, there is a 'Register' button, a radio button for 'I have read and agree to the Privacy Policy', and a link for 'I already have an account'.

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



1. Front thumb electrodes
2. Handle with retractable cable
3. Colour LCD scrolling display
4. Plastic base
5. Tempered glass tray
6. Stainless steel foot electrodes
7. Rear electrodes for other fingers
8. Silicone feet
9. Start button
10. USB-C socket for charging

To download the Fitdays+ app, scan the QR code above
 Mobile phones supported: from **iOS 8.0** or **Android 6.0**

1. PRECAUTIONS FOR USE

- ▶ Do not use when wearing implants or a pacemaker.
- ▶ Not for use by pregnant women
- ▶ The measurements indicated are to be considered as indicators but do not replace medical advice
- ▶ Consult your doctor before any major change in diet and/or physical activity
- ▶ Only step on the scales with feet that have been thoroughly dried (risk of slipping)
- ▶ Do not climb on the edge of the scale, do not jump on the scale
- ▶ Do not overload (180 kg max)
- ▶ Do not shock the appliance or drop an object onto it
- ▶ Do not immerse in liquid
- ▶ Do not use detergent for cleaning
- ▶ Place the appliance on a hard, flat, dry surface
- ▶ Measurements unreliable or impossible < 18 years (imperfect bone growth)

2. FIRST STEPS

Remove the transport lock (prevents unloading during transport and storage) by pressing and holding the button (1) - also used to reboot.
 Make sure the device is sufficiently charged; if necessary, charge it by plugging it into the USB-C socket (2).

To enable the device to calculate your values, you must first enter your personal details in the application. There are then 2 ways of using the device:

- a. Quick-mode: without using the application, your main data scrolls across the device screen
- b. Connected mode: once your phone is connected, the application displays all your body analysis data.

NB: the 'Account' section on the home page allows you to synchronise with third-party applications such as Apple Health, Fitbit, etc.


The automatic recognition function automatically identifies users with a weight difference of more than 2 kg.

The symbols on the scale display are as follows:

1. Battery level: flashes when you need to recharge
2. Bluetooth connection
3. Large digits display measured values

The scrolling dotted lines indicate that the weighing process is in progress, then 3 x flashes indicate that the weight value is locked/saved.

The device then switches to scrolling screen mode:

3. VOCAL SUPPORT

- ▶ Available in English and French
- ▶ Weight measurement
- ▶ Fat measurement
- ▶ Correct handlebar position and grip

4. USING THE FITDAYS APPLICATION

- ▶ Download the application (see QR-code on front page)
- ▶ Activate the phone's bluetooth and wait for the connection, the scale will automatically enter search mode
- ▶ Once the connection has been established, the corresponding icon will appear on the scale screen
- ▶ In the application, choose the scale / My scale (you can then give it a name of your choice)
- ▶ Register with the application, then enter your personal details so that the application can perform its calculations

A screenshot of the 'Register' screen in the Fitdays application. At the top, it shows the time '08:42' and signal/battery icons. Below the title 'Register' is a 'Feedback' link. The form has two input fields: 'Email' with a placeholder 'Please enter your Email' and 'Password' with a placeholder 'Please input password' and a strength indicator. Below the password field is a list of requirements: 'password should contain 8-16 characters' and 'Include at least two kinds among numbers, uppercase and lowercase letters, symbols'. At the bottom is a 'Register' button, a radio button for 'I have read and agree to the Privacy Policy', and a link for 'I already have an account'.

Press down on the scale with your foot to activate it. Wait until the display shows '0.0' kg.

Start the application on your phone, then step onto the scale barefoot while holding the handle correctly, with your arms at a 45° angle to your body.

In just a few seconds, you'll be able to see a complete analysis of your body values, position them in relation to average values, track their evolution and draw up a detailed report that you can print out directly. The QR-code below gives you a demonstration of the main uses:



You can add users by going to the menu Device / user management / add user

You can also share the device remotely by going to the menu Device / share device / share

5. ERRORS / FAQ

- ▶ Hold the arms at an approximate angle of 45° to the body.
- ▶ The fingers must be correctly aligned on the electrodes with which they must all be in contact

- ▶ The application must be connected and the personal data filled in, otherwise the device will only be able to give the weight
- ▶ It is essential to be barefoot (no socks) and dry
- ▶ Do not descend too quickly so that the device can complete its analysis cycle (10 sec. after the weight has stopped flashing)
- ▶ Be careful not to mix users in the application at the time of body analysis
- ▶ Incorrect weight: check that you are on hard, flat ground (possibly do a reboot with button 9)
- ▶ Incorrect ratios: check the accuracy of the parameters you have entered in the app.
- ▶ No bluetooth connection:
 - o the scale must be within 10 m. of the phone.
 - o it may have connected to another device by accident
 - o the phone must be running at least iOS 8.0 or Android 6.0
 - o if necessary, press off/on on your phone's bluetooth function