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
Bodyfat Scale LB8838 - Instruction manual

LB8838IMINTERV1

Thank you for placing your trust in Little Balance.


Although our device is very intuitive to use, we have made every effort to provide clear and comprehensive instructions.

IMPORTANT: Before using the device, please read these instructions carefully. Failure to comply with these instructions may result in personal injury or damage to the unit.

To download the Fitdays+ app, scan the QR-code
Cell phones supported: from iOS 8.0 or Android 6.0

You can contact Littlebalance at balancy@littlebalance.fr
The development team can be reached via the app: Account / feedback / scale.

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PRECAUTIONS FOR USE

- ▶ Do not use when wearing implants or a pacemaker
- ▶ Not for use by pregnant women
- ▶ The measurements indicated are to be considered as indicators but do not replace medical advice
- ▶ Consult your doctor before any major change in diet and/or physical activity
- ▶ Stand on the scale only with feet thoroughly dried (risk of slipping)
- ▶ Do not climb on the edge of the scale, do not jump on the scale
- ▶ Do not overload (180 kg max)
- ▶ Do not shock the device or drop an object onto it
- ▶ Do not immerse in liquid
- ▶ Do not use detergents for cleaning
- ▶ Place the device on a hard, flat, dry surface
- ▶ Measurements unreliable or impossible < 18 years (imperfect bone growth)

FEATURES OF YOUR DEVICE

Your device has several features that make it very practical to use:

1. It has a double display:

- the upper screen displays your weight
- the lower screen displays your body analysis:
- 8 analyses directly readable from the scale display

2. It connects to your phone via Bluetooth, not Wi-Fi:

simple and practical, this feature eliminates the need to connect via the wifi network and its access codes.

3. Once it has been connected to the application for the first time and you have entered your personal data, it no longer needs to be systematically connected to carry out your body analysis because no longer needs to be systematically logged in to perform your body analysis, as:

- it has an automatic offline identification system that enables it to recognize users (provided they have a minimum weight difference of 2 kg between them)
- it has a memory function that allows 24 weighings per user to be recorded and phone when the application and connection are activated
- So, there's no need to carry your phone with you every time you do a body analysis:
 - 8 parameters can be read directly and immediately on the scale
 - data are memorized and sent to the application the next time you connect.

FIRST STEPS

Remove the transport lock (prevents unloading during transport and storage) by pressing and holding the button (1) - also used for rebooting

Make sure the device is sufficiently charged; if necessary, charge it by plugging into the USB-C socket (2)

To enable the device to calculate your values, you must first enter your personal data in the application. There are then 2 ways of using the device:

- a. Quick mode : without using the application, your main data
 - i. 8 data items are displayed on the device screen for immediate readability without having your phone in hand
 - ii. all data are stored in memory, to be sent to the application the next time you log on.
- b. Connected mode: once your phone is connected, the application displays all your body analysis data.

NB: the "Account" section on the home page enables synchronization with third-party applications such as Apple Health, Fitbit...

The automatic recognition function automatically identifies users with a weight difference of over 2 kg.

The scale displays are as follows:

1. top screen
 - a. bluetooth connection symbols
 - b. weight with a precision of 50 gr.
2. lower screen
 - a. heart rate
 - b. body fat
 - c. body mass index (BMI)
 - d. water content
 - e. bone mass
 - f. muscle mass
 - g. body type

And, of course, the Bluetooth connection and low battery symbols!


FOR SMOOTH OPERATION

The application must be connected and personal data filled in, otherwise the device will only give the weight

- ▶ You must be barefoot (no socks) and dry.
- ▶ Do not descend too quickly for the device to complete its analysis cycle (10 sec. after the weight has stopped flashing)
- ▶ Be careful not to mix users in the application at the time of body analysis
- ▶ Incorrect weight: check that you're on hard, flat ground (possibly do a reboot with button 9)
- ▶ Incorrect ratios: check the accuracy of the parameters you've entered in the app
- ▶ No Bluetooth connection:
 - the scale must be within 10 m. of the phone
 - it may have connected to another device by accident
 - the phone must be running at least iOS 8.0 or Android 6.0
 - if necessary, press off/on on your phone's bluetooth function

USING THE FITDAYS APPLICATION

- ▶ Download the app (see QR-code on front page)
- ▶ Activate the phone's Bluetooth and wait for the connection, the scale automatically enters search mode
- ▶ Once the connection is established, the corresponding icon appears on the scale display
- ▶ In the application, select the scale / My scale (you can then give it a name of your choice)
- ▶ Register with the application, then enter your personal data so that the application can perform its calculations

Press down on the scale with your foot to activate it. Wait until the display shows “0.0” kg.

Start the app on your phone, then step onto the scale barefoot.

In just a few seconds, the device performs a complete analysis of your body values, and the app positions them in relation to average values. You can follow the evolution of these values, data by data, and even draw up a detailed report that you can print out directly. The QR-code below gives you a demonstration of the main uses:

You can add users by going to the menu Device / user management / add user

You can also share the device remotely by going to menu Device / share device / share

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