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
Bodyfat scale LB8861 - Instructions for use

LB8861IMINTERV1

Thank you for placing your trust in Little Balance.

Although our equipment is very intuitive to use, we have endeavoured to provide clear and comprehensive instructions.

IMPORTANT: Before using the appliance, please read the instructions in this manual carefully. Failure to comply with the instructions or this manual may result in personal injury or damage to the appliance.


To download the Fitdays+ app, scan the QR code

Mobile phones supported: from iOS 8.0 or Android 6.0

You can contact Littlebalance at balancy@littlebalance.fr

The development team can be reached via the app: Account / feedback / scales

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1. PRECAUTIONS FOR USE

- ▶ Do not use if you have implants or a pacemaker.
- ▶ Do not use by pregnant women
- ▶ The measurements indicated are to be considered as indicators but do not replace medical advice
- ▶ Consult your doctor before any major change in diet and/or physical activity
- ▶ Only step on the scales with feet that have been thoroughly dried (risk of slipping)
- ▶ Do not stand on the edge of the scale, do not jump on the scale
- ▶ Do not overload (180 kg max)
- ▶ Do not shock the appliance or drop an object onto it
- ▶ Do not immerse in liquid
- ▶ Do not use detergent for cleaning
- ▶ Place the appliance on a hard, flat, dry surface
- ▶ Measurements unreliable or impossible < 18 years (imperfect bone growth)

2. FEATURES OF YOUR APPLIANCE

Your appliance has a number of features that make it very practical to use:

1. It has a central indicator light that displays

colour-coded changes in weight:

- blue = stable weight
- green = weight loss
- red = weight gain

Warning: when the device is not connected, the user is only identified is only performed if the weight difference is < 2kg.

If the difference is > 2kg, automatic identification will not take place and the default colour code will come on.

2. It connects to your phone via Bluetooth rather than Wi-Fi:

simple and practical, this feature avoids having to connect via the wifi network and its access codes.

3. Once it has been connected to the application for the first time and you have entered your personal data, it will no longer need to be systematically logged in to carry out your body analyses because :

- a. it has an automatic offline identification system which enables it to recognise users (provided they are at least 2 kg different in weight)
- b. it has a memory function that allows it to record 24 weighings per user, and then send them to the phone when the application and connection are activated
- c. So there's no need to carry your phone with you every time you do a body analysis:
- e. 7 parameters can be read directly and immediately on the scale
- f. the data is stored and sent to the application the next time you connect

3. FIRST STEPS

Remove the transport lock (prevents unloading during transport and storage) by pressing and holding the button (1) - also used to reboot. Make sure the device is sufficiently charged; if necessary, charge it by plugging it into the USB-C socket (2).

To enable the device to calculate your values, you must first enter your personal details in the application. There are then 2 ways of using the device:

- Quick-mode: without using the application, your main data
- 7 data items are displayed on the device's screen for immediate readability without having your phone in your hand.
- All data is stored in memory and sent to the application the next time you connect
- Connected mode: once your phone is connected, the application displays all your body analyses on it

NB: the 'Account' section on the home page allows you to synchronise with third-party applications such as Apple Health, Fitbit, etc.

The automatic recognition function automatically identifies users with a weight difference of more than 2 kg.

The scale displays are as follows:

1. Weight: indicated with an accuracy of 50 g.
2. Bluetooth connection icon
3. Low battery icon
4. Illuminated central display showing weight trend
 - blue = stable weight
 - green = weight loss
 - red = weight gain


4. FOR CORRECT OPERATION

- ▶ The application must be connected and the personal data filled in, otherwise the device will only be able to give the weight
- ▶ It is essential to be barefoot (no socks) and dry
- ▶ Do not descend too quickly so that the device can complete its analysis cycle (10 sec. after the weight has stopped flashing)
- ▶ Be careful not to mix users in the application at the time of body analysis
- ▶ Incorrect weight: check that you are on hard, flat ground (possibly do a reboot with button 9)
- ▶ Incorrect ratios: check the accuracy of the parameters you have entered in the app.
- ▶ No bluetooth connection:
 - the scale must be within 10 m. of the phone.
 - it may have connected to another device by accident
 - the phone must be running at least iOS 8.0 or Android 6.0
 - if necessary, press off/on on your phone's Bluetooth function

5. USING THE FITDAYS APPLICATION

- ▶ Download the application (see QR-code on front page)
- ▶ Activate the phone's bluetooth and wait for the connection, the scale will automatically enter search mode
- ▶ Once the connection has been established, the corresponding icon will appear on the scale screen
- ▶ In the application, choose the scale / My scale (you can then give it a name of your choice)
- ▶ Register with the application, then enter your personal details so that the application can perform its calculations

Press down on the scale with your foot to activate it. Wait until the display shows 0.0 kg.

Start the app on your phone, then step onto the scale barefoot.

In just a few seconds, the device will carry out a complete analysis of your body values and the application will position them in relation to average values. You can follow the evolution of your values, data by data, and even draw up a detailed report that you can print out directly. The QR-code below gives you a demonstration of the main uses:

You can add users by going to the menu Device / user management / add user

You can also share the device remotely by going to the menu Device / share device / share

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